## azœka

## **MURAPR BULK** Sample Collection & Storage Media (RUO)



## **PRODUCT DESCRIPTION**

Conveniently collect, store and transport your samples with mWRAPR BULK Sample collection & storage media. mWRAPR BULK sample collection & media storage potently inactivates pathogens and arrests their growth while effectively preserving genetic material for downstream processing.

## **FEATURES**

- No liquid nitrogen cold storage required
  - Protection from freeze-thaw cycles
  - Long-term storage at -20°C
  - optional) Storage at -80°C
- Easy sample processing
  - DNA / RNA can be directly isolated from mWRAPR media using common commercial extraction kits
- Stabilize genetic material for up to 14 days at room temperature
- Zero contamination shelf-life for up to 2 years
- Accurate and reproducible DNA/RNA yields for downstream processing & analysis

Product Name	Cat. No.	Sample Input Vol.	Stabilization solution Vol.	Pack Size	Price € (Euros)
mWRAPR BULK					